

## Urine Drug Screening Cutoffs

Options Lab, Inc. utilizes the following drug screening methods: EIA (Enzyme Immunoassay) and LC/MS (Liquid Chromatography/Mass Spectrometry)

DRUG CLASS	CUT-OFF
6-Acetylmorphine (6-AM)	10 ng/ml
Amphetamines	500 ng/ml
Barbiturates	200 ng/ml
Benzodiazepines	200 ng/ml
Buprenorphine	5 ng/ml
Cannabinoids (Marijuana)	20 or 50 ng/ml
Cocaine	150 ng/ml
Cotinine	250 ng/ml
Dextromethorphan	100 ng/ml
Ethanol	20 mg/dL
Ethyl Glucuronide (EtG)	100 or 500 ng/ml
Ethyl Sulfate (EtS)	25 or 125 ng/ml
Fentanyl	1 ng/ml
Methadone	300 ng/ml
Opiates	300 ng/ml
Oxycodone	100 ng/ml
Phencyclidine (PCP)	25 ng/ml
Synthetic Cannabinoids (K2-1)	10 ng/ml
Synthetic Cannabinoids (K2-2)	10 ng/ml